

PRODUCT LINE-UP

Safety Made Simple for Every Setting

From first aid and professional rescue to everyday safety and sports recovery, Keepu makes safety and recovery simpler for every setting.



01 FACILITY SAFETY

Twist Care AS-Basic200

A reusable air splint designed for fast, intuitive use by non-professionals.

APPROX. 10 SEC | 4 BODY AREAS | CO₂ INFLATION



02 PROFESSIONAL RESCUE

Airfix AS-Pro200

A reusable air splint optimized for repeated use in professional rescue settings.

AUTO INFLATION | 300 CYCLES | REUSABLE



03 EVERYDAY SAFETY

Tight Care AS-Lite100

An intuitive air splint designed to make everyday safety easier and more accessible.

WITHIN 20 SEC | NO EQUIPMENT | REUSABLE



04 SPORTS RECOVERY

ARI

A sports recovery product that supports during activity and cools right after.

KNEE SUPPORT | NO FREEZER | 20-MIN COOLING

Three Air Splints for Different Settings

Built on the same pneumatic stabilization principle,
each model is optimized for facility use,
professional rescue, or everyday safety.



Twist Care AS-Basic200

FACILITY SAFETY

CO₂ cartridge inflation
for rapid first response



Airfix AS-Pro200

PROFESSIONAL RESCUE

Automatic open-type design
for repeated professional use



Tight Care AS-Lite100

EVERYDAY SAFETY

Simple, equipment-free design
for everyday preparedness

Category	Twist Care	Airfix	Tight Care
Primary Use	Facility Readiness	Professional Rescue	Everyday Safety Preparedness
Inflation Method	CO ₂ Cartridge	Automatic Air Pump	Integrated Oral Inflation Port
Inflation Time	Approx. 10 sec	Approx. 10 sec	Within 20 sec
Recommended Pressure	Approx. 25–40 mmHg	Approx. 25–40 mmHg	Approx. 25–35 mmHg
Reusable	Yes	Yes	Yes
Key Feature	Fast, Intuitive Response	Easy Cleaning & Reuse	No Extra Equipment
Regulatory Status	Reusable Medical Device — Registered	Reusable Medical Device — Registered	Reusable Medical Device — Registration Planned

MATERIAL TECHNOLOGY

Air Splint Material Co-Developed with LG Chem

Co-developed with LG Chem to withstand repeated inflation and real-world use, with material properties optimized for pneumatic air splints.

FLEXIBILITY

Flexible Fit

DURABILITY

Designed for Repeated Use

MATTE SURFACE

Enhanced Comfort

4H

40 mmHg Pressure Retention

300×

Airfix Inflation / Deflation Test

200mmHg

Structural Safety Test

PATENT 2


2 Patent Applications Filed

Support While You Move. Ice Right After.

Wear it as a knee support during activity, then activate the instant ice pack to cool the target area immediately after exercise.




ARI — Knee Support + Instant Ice Packs



01

Knee Support During Activity

Elastic support provides a secure, comfortable fit during movement.



02

Instant Cooling, No Pre-Freezing

Activate the instant ice pack and use it immediately after exercise.



03

Approx. 20-Min Targeted Cooling

Insert the ice pack into the pocket to cool the target area hands-free.

20 MIN

Cooling Duration

NO FREEZER

No Pre-Freezing Required

M / L

Two Sizes

REUSABLE

Washable & Reusable Support

PACKAGE CONTENTS

Knee Support
x 1

Instant Ice Packs
x 10

Knee Support

Washable & Reusable

Ice Pack

Single Use · Refills Available

Product Category

General Sports Product

SIZE GUIDE

SIZE	TOP WIDTH	BOTTOM WIDTH	LENGTH
M	17 cm	15 cm	27 cm
L	18.5 cm	16.5 cm	27 cm

HOMEPAGE
keepu.pe.kr

E-MAIL
contact@team-keepu.com

TEL
(+82)070-8065-1291

ADDRESS
55F, Busan International Finance Center
40 Munhyeongeumyung-ro, Nam-gu
Busan, Republic of Korea

